


## TC HEANEY NEEDLE HOLDER

curved; serrated

26.0cm (10.4")



**SKU:** 051-141-26

**Categories:** [Open Surgery](#), [Suture](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>KLS Martin</b>	20-643-26
<b>AUTHORS</b>	<a href="#">HEANEY</a>
<b>SPECIFICATIONS</b>	<a href="#">curved</a> , <a href="#">serrated</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>KLS Martin</b>	20-643-26
<b>AUTHORS</b>	<a href="#">HEANEY</a>

## SPECIFICATIONS

[curved](#), [serrated](#)