

TC HEGAR NEEDLE HOLDER

straight; serrated

20.0cm (8")

SKU: 051-084-20

Categories: [Open Surgery](#), [Suture](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM077
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-662-20
AUTHORS	HEGAR
Medicon	10.18.71
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM077

MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-662-20
AUTHORS	HEGAR
Medicon	10.18.71
SPECIFICATIONS	serrated , straight