


TC HEGAR-VASCULAR NEEDLE HOLDER

straight; serrated

18.0cm (7.2")



SKU: 051-090-18

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
AUTHORS	HEGAR, VASCULAR
SPECIFICATIONS	serrated, straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
AUTHORS	HEGAR, VASCULAR
SPECIFICATIONS	serrated, straight