


## TC HEGAR-VASCULAR NEEDLE HOLDER

straight; serrated

20.0cm (8")



**SKU:** 051-090-20

**Categories:** [Open Surgery](#), [Suture](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>AUTHORS</b>	<a href="#">HEGAR, VASCULAR</a>
<b>Medicon</b>	10.25.71
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>AUTHORS</b>	<a href="#">HEGAR, VASCULAR</a>
<b>Medicon</b>	10.25.71

## SPECIFICATIONS

[serrated](#), [straight](#)