


TC HEGAR-VASCULAR NEEDLE HOLDER

straight; serrated

14.0cm (5.6")



SKU: 051-090-14

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Needle Holders |
| MATERIAL | Tungsten Carbide(TC) |
| AUTHORS | HEGAR, VASCULAR |
| Medicon | 10.25.65 |
| SPECIFICATIONS | serrated , straight |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------|--------------------------------------|
| UNIT | pcs. |
| TYPE | Needle Holders |
| MATERIAL | Tungsten Carbide(TC) |
| AUTHORS | HEGAR, VASCULAR |
| Medicon | 10.25.65 |

SPECIFICATIONS

[serrated](#), [straight](#)