

TC LANGENBECK NEEDLE HOLDER

straight; serrated; without catch

20.0cm (8")

SKU: 051-360-20

Categories: [Open Surgery](#), [Suture](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM193
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-694-20
AUTHORS	LANGENBECK
Medicon	11.36.70
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM193

MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-694-20
AUTHORS	LANGENBECK
Medicon	11.36.70
SPECIFICATIONS	serrated , straight