

TC LANGENBECK NEEDLE HOLDER

straight; serrated; without catch

23.0cm (9.2")

SKU: 051-360-23

Categories: [Open Surgery](#), [Suture](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
AUTHORS	LANGENBECK
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
AUTHORS	LANGENBECK
SPECIFICATIONS	serrated , straight