

TC LICHTENBERG NEEDLE HOLDER

straight; serrated; with catch

17.0cm (6.8")

SKU: 051-366-17

Categories: [Open Surgery](#), [Suture](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
AUTHORS	LICHTENBERG
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
AUTHORS	LICHTENBERG
SPECIFICATIONS	serrated , straight