


TC MARTINI NEEDLE HOLDER

straight; smooth

14.0cm (5.6")


SKU: 051-120-14

Categories: [Open Surgery](#), [Suture](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
AUTHORS	MARTINI
SPECIFICATIONS	straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
AUTHORS	MARTINI
SPECIFICATIONS	straight