

## TC MATHIEU-KOCHER NEEDLE HOLDER

straight; serrated; with catch

17.0cm (6.8")

**SKU:** 051-358-17

**Categories:** [Open Surgery](#), [Suture](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>Aesculap</b>	BM150
<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>KLS Martin</b>	20-674-17
<b>AUTHORS</b>	<a href="#">KOCHER</a> , <a href="#">MATHIEU</a>
<b>Medicon</b>	11.18.67
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>Aesculap</b>	BM150

<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>KLS Martin</b>	20-674-17
<b>AUTHORS</b>	<a href="#">KOCHER, MATHIEU</a>
<b>Medicon</b>	11.18.67
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>