

TC MATHIEU-KOCHER NEEDLE HOLDER

straight; serrated; with catch

17.0cm (6.8")



SKU: 051-358-17

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Needle Holders |
| Aesculap | BM150 |
| MATERIAL | Tungsten Carbide(TC) |
| KLS Martin | 20-674-17 |
| AUTHORS | KOCHER , MATHIEU |
| Medicon | 11.18.67 |
| SPECIFICATIONS | serrated , straight |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------|--------------------------------|
| UNIT | pcs. |
| TYPE | Needle Holders |
| Aesculap | BM150 |

| | |
|-----------------------|---|
| MATERIAL | Tungsten Carbide(TC) |
| KLS Martin | 20-674-17 |
| AUTHORS | KOCHER, MATHIEU |
| Medicon | 11.18.67 |
| SPECIFICATIONS | serrated , straight |