

## TC MATHIEU NEEDLE HOLDER

straight; serrated; with catch

17.0cm (6.8")

**SKU:** 051-350-17

**Categories:** [Open Surgery](#), [Suture](#)



### ZUSÄTZLICHE INFORMATIONEN

|                       |   |
|-----------------------|---|
| <b>UNIT</b>           | <a href="#">pcs.</a>                                |
| <b>TYPE</b>           | <a href="#">Needle Holders</a>                      |
| <b>Aesculap</b>       | BM366   |
| <b>MATERIAL</b>       | <a href="#">Tungsten Carbide(TC)</a>                |
| <b>KLS Martin</b>     | 20-676-17   |
| <b>AUTHORS</b>        | <a href="#">MATHIEU</a>                             |
| <b>Medicon</b>        | 11.20.67  |
| <b>SPECIFICATIONS</b> | <a href="#">serrated</a> , <a href="#">straight</a> |



### ZUSÄTZLICHE INFORMATIONEN

|                 |                                |
|-----------------|--------------------------------|
| <b>UNIT</b>     | <a href="#">pcs.</a>           |
| <b>TYPE</b>     | <a href="#">Needle Holders</a> |
| <b>Aesculap</b> | BM366                          |

|                       |   |
|-----------------------|---|
| <b>MATERIAL</b>       | <a href="#">Tungsten Carbide(TC)</a>                |
| <b>KLS Martin</b>     | 20-676-17   |
| <b>AUTHORS</b>        | <a href="#">MATHIEU</a>                             |
| <b>Medicon</b>        | 11.20.67  |
| <b>SPECIFICATIONS</b> | <a href="#">serrated</a> , <a href="#">straight</a> |