

TC MATHIEU NEEDLE HOLDER

straight; serrated; with catch

20.0cm (8")

SKU: 051-350-20

Categories: [Open Surgery](#), [Suture](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM154
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-676-20
AUTHORS	MATHIEU
Medicon	11.20.70
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM154

MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-676-20
AUTHORS	MATHIEU
Medicon	11.20.70
SPECIFICATIONS	serrated , straight