

TC MAYO-HAUBERISSER WIRE TWISTER

straight; serrated

15.0cm (6")



SKU: 051-244-15

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Wire twisting Forceps
MATERIAL	Tungsten Carbide(TC)
AUTHORS	HAUBERISSER , MAYO
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Wire twisting Forceps
MATERIAL	Tungsten Carbide(TC)
AUTHORS	HAUBERISSER , MAYO
SPECIFICATIONS	serrated , straight