

## TC MAYO-HAUBERISSER WIRE TWISTER

straight; serrated

15.0cm (6")



**SKU:** 051-244-15

**Categories:** [Open Surgery](#), [Suture](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Wire twisting Forceps</a>
<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>AUTHORS</b>	<a href="#">HAUBERISSER</a> , <a href="#">MAYO</a>
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Wire twisting Forceps</a>
<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>AUTHORS</b>	<a href="#">HAUBERISSER</a> , <a href="#">MAYO</a>
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>