


TC MAYO-HEGAR NEEDLE HOLDER

straight; serrated

18.0cm (7.2")



SKU: 051-080-18

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM066
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-642-18
AUTHORS	HEGAR , MAYO
Medicon	10.18.68
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM066

MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-642-18
AUTHORS	HEGAR, MAYO
Medicon	10.18.68
SPECIFICATIONS	serrated, straight