


## TC MAYO-HEGAR NEEDLE HOLDER

straight; serrated

24.0cm (9.6")




**SKU:** 051-080-24

**Categories:** [Open Surgery](#), [Suture](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>Aesculap</b>	BM078
<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>KLS Martin</b>	20-642-24
<b>AUTHORS</b>	<a href="#">HEGAR</a> , <a href="#">MAYO</a>
<b>Medicon</b>	10.18.74
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>Aesculap</b>	BM078

<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>KLS Martin</b>	20-642-24
<b>AUTHORS</b>	<a href="#">HEGAR, MAYO</a>
<b>Medicon</b>	10.18.74
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>