

## TC MAYO-HEGAR NEEDLE HOLDER LH

straight; serrated; for left hand use

18.0cm (7")

**SKU:** 051-080-18/LH

**Categories:** [Open Surgery](#), [Suture](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>AUTHORS</b>	<a href="#">HEGAR</a> , <a href="#">MAYO</a>
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>AUTHORS</b>	<a href="#">HEGAR</a> , <a href="#">MAYO</a>
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>