


## TC MICRO RYDER NEEDLE HOLDER

straight; serrated

15.0cm (6")



**SKU:** 051-132-15

**Categories:** [Open Surgery](#), [Suture](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>Aesculap</b>	BM055
<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>KLS Martin</b>	20-654-15
<b>AUTHORS</b>	<a href="#">MICRO</a> , <a href="#">RYDER</a>
<b>Medicon</b>	10.28.65
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>Aesculap</b>	BM055

<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>KLS Martin</b>	20-654-15
<b>AUTHORS</b>	<a href="#">MICRO</a> , <a href="#">RYDER</a>
<b>Medicon</b>	10.28.65
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>