


TC MICRO RYDER NEEDLE HOLDER

straight; serrated

13.0cm (5.2")



SKU: 051-132-13

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM054
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-654-13
AUTHORS	MICRO , RYDER
Medicon	10.28.63
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM054

MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-654-13
AUTHORS	MICRO , RYDER
Medicon	10.28.63
SPECIFICATIONS	serrated , straight