


TC MICROVASCULAR NEEDLE HOLDER

straight; serrated

18.0cm (7.2")



SKU: 051-122-18

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-666-02
AUTHORS	MICROVASCULAR
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-666-02
AUTHORS	MICROVASCULAR

SPECIFICATIONS

[serrated](#), [straight](#)