


TC MICROVASCULAR NEEDLE HOLDER

straight; serrated

20.0cm (8")



SKU: 051-122-20

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-666-03
AUTHORS	MICROVASCULAR
Medicon	10.24.17
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-666-03

AUTHORS

[MICROVASCULAR](#)

Medicon

10.24.17

SPECIFICATIONS

[serrated](#), [straight](#)