


## TC MICROVASCULAR NEEDLE HOLDER

straight; serrated

20.0cm (8")



**SKU:** 051-122-20

**Categories:** [Open Surgery](#), [Suture](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>KLS Martin</b>	20-666-03
<b>AUTHORS</b>	<a href="#">MICROVASCULAR</a>
<b>Medicon</b>	10.24.17
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>KLS Martin</b>	20-666-03

**AUTHORS**

[MICROVASCULAR](#)

**Medicon**

10.24.17

**SPECIFICATIONS**

[serrated](#), [straight](#)