


TC MICROVASCULAR NEEDLE HOLDER

straight; serrated

14.0cm (5.6")



SKU: 051-122-14

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-666-00
AUTHORS	MICROVASCULAR
Medicon	10.24.65
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-666-00

AUTHORS

[MICROVASCULAR](#)

Medicon

10.24.65

SPECIFICATIONS

[serrated](#), [straight](#)