

## TC NEEDLE HOLDER Ø 3MM, CURVED RIGHT

with Aluminium handle, teflon coated

WL 290mm, LL-Connection



**SKU:** 614-115-33

**Categories:** [Endoscopy](#), [Laparoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Needle Holders</a> , <a href="#">Pediatric Laparoscopy</a>
<b>TYPE</b>	<a href="#">LOFTLINE</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Needle Holders</a> , <a href="#">Pediatric Laparoscopy</a>
<b>TYPE</b>	<a href="#">LOFTLINE</a>