

## TC NEEDLE HOLDER Ø 3MM, STRAIGHT

with Aluminium handle, teflon coated

WL 290mm, LL-Connection



**SKU:** 614-100-33

**Categories:** [Endoscopy](#), [Laparoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

|              |  |
|--------------|--|
| <b>UNIT</b>  | <a href="#">pcs.</a>   |
| <b>RANGE</b> | <a href="#">Needle Holders</a> , <a href="#">Pediatric Laparoscopy</a> |
| <b>TYPE</b>  | <a href="#">LOFTLINE</a>   |

### ZUSÄTZLICHE INFORMATIONEN

|              |  |
|--------------|--|
| <b>UNIT</b>  | <a href="#">pcs.</a>   |
| <b>RANGE</b> | <a href="#">Needle Holders</a> , <a href="#">Pediatric Laparoscopy</a> |
| <b>TYPE</b>  | <a href="#">LOFTLINE</a>   |