

## TC NEEDLE HOLDER Ø 5MM, CURVED RIGHT

with Aluminium handle, teflon coated

WL 330mm, LL-Connection



**SKU:** 614-115-53

**Categories:** [Endoscopy](#), [Laparoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Needle Holders](#)

**TYPE** [LOFTLINE](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Needle Holders](#)

**TYPE** [LOFTLINE](#)