

TC NEEDLE HOLDER Ø 5MM, CURVED RIGHT

with Aluminium handle, teflon coated

WL 330mm, LL-Connection



SKU: 614-115-53

Categories: [Endoscopy](#), [Laparoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Needle Holders](#)

TYPE [LOFTLINE](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Needle Holders](#)

TYPE [LOFTLINE](#)