

## TC NEEDLE HOLDER; STRAIGHT

axial handle straight; with ratchet

Luer-Lock; Ø 5.0mm; 330mm



**SKU:** 614-405-53

**Categories:** [Endoscopy](#), [Laparoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Needle Holders](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Needle Holders](#)