


TC OLSEN-HEGAR NEEDLE HOLDER

straight; serrated; with scissors

17.0cm (6.8")



SKU: 051-310-17

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM124
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-668-17
AUTHORS	HEGAR , OLSEN
Medicon	10.60.67
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM124

MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-668-17
AUTHORS	HEGAR, OLSEN
Medicon	10.60.67
SPECIFICATIONS	serrated , straight