

## TC OLSEN-HEGAR NEEDLE HOLDER

straight; serrated; with scissors

19.0cm (7.6")



**SKU:** 051-310-19

**Categories:** [Open Surgery](#), [Suture](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>KLS Martin</b>	20-672-18
<b>AUTHORS</b>	<a href="#">HEGAR</a> , <a href="#">OLSEN</a>
<b>Medicon</b>	10.60.69
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>KLS Martin</b>	20-672-18

**AUTHORS**

[HEGAR](#), [OLSEN](#)

**Medicon**

10.60.69

**SPECIFICATIONS**

[serrated](#), [straight](#)