

TC OLSEN-HEGAR NEEDLE HOLDER

straight; serrated; with scissors

14.0cm (5.6")




SKU: 051-310-14

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Needle Holders |
| Aesculap | BM128 |
| MATERIAL | Tungsten Carbide(TC) |
| KLS Martin | 20-672-14 |
| AUTHORS | HEGAR , OLSEN |
| Medicon | 10.60.64 |
| SPECIFICATIONS | serrated , straight |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------|--------------------------------|
| UNIT | pcs. |
| TYPE | Needle Holders |
| Aesculap | BM128 |

| | |
|-----------------------|---|
| MATERIAL | Tungsten Carbide(TC) |
| KLS Martin | 20-672-14 |
| AUTHORS | HEGAR, OLSEN |
| Medicon | 10.60.64 |
| SPECIFICATIONS | serrated , straight |