

TC POTTS-SMITH DRESSING FORCEPS

straight; serrated

18.0cm (7.2")

SKU: 021-026-18
 Categories: [Open Surgery](#), [Forceps](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD156
MATERIAL	Tungsten Carbide(TC)
KLS Martin	12-258-18
AUTHORS	POTTS , SMITH
Medicon	06.02.68
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD156

MATERIAL	Tungsten Carbide(TC)
KLS Martin	12-258-18
AUTHORS	POTTS , SMITH
Medicon	06.02.68
SPECIFICATIONS	serrated , straight