

TC POTTS-SMITH DRESSING FORCEPS

straight; serrated

15.0cm (6")

SKU: 021-026-15
 Categories: [Open Surgery](#), [Forceps](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD154
MATERIAL	Tungsten Carbide(TC)
KLS Martin	12-258-15
AUTHORS	POTTS , SMITH
Medicon	06.02.65
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD154

MATERIAL	Tungsten Carbide(TC)
KLS Martin	12-258-15
AUTHORS	POTTS , SMITH
Medicon	06.02.65
SPECIFICATIONS	serrated , straight