


TC RYDER "FRENCH EYE" NEEDLE HOLDER

straight; serrated

23.0cm (9.2")



SKU: 051-134-23

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-650-23
AUTHORS	RYDER
Medicon	10.26.73
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-650-23

AUTHORS

[RYDER](#)

Medicon

10.26.73

SPECIFICATIONS

[serrated](#), [straight](#)