


TC SAROT NEEDLE HOLDER

straight; serrated

20.0cm (8")



SKU: 051-224-20

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Needle Holders |
| MATERIAL | Tungsten Carbide(TC) |
| AUTHORS | SAROT |
| Medicon | 10.22.72 |
| SPECIFICATIONS | serrated , straight |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------|--------------------------------------|
| UNIT | pcs. |
| TYPE | Needle Holders |
| MATERIAL | Tungsten Carbide(TC) |
| AUTHORS | SAROT |
| Medicon | 10.22.72 |

SPECIFICATIONS

[serrated](#), [straight](#)