


TC SAROT NEEDLE HOLDER

straight; serrated

26.0cm (10.4")



SKU: 051-224-26

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM023
MATERIAL	Tungsten Carbide(TC)
KLS Martin	26-640-26
AUTHORS	SAROT
Medicon	10.22.77
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM023

MATERIAL	Tungsten Carbide(TC)
KLS Martin	26-640-26
AUTHORS	SAROT
Medicon	10.22.77
SPECIFICATIONS	serrated , straight