


TC SAROT NEEDLE HOLDER

straight; serrated

18.0cm (7.2")




SKU: 051-224-18

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM020
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-640-18
AUTHORS	SAROT
Medicon	10.22.68
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM020

MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-640-18
AUTHORS	SAROT
Medicon	10.22.68
SPECIFICATIONS	serrated , straight