

TC SEMKEN DRESSING FORCEPS

straight; serrated

18.0cm (7.2")


SKU: 021-040-18

Categories: [Open Surgery](#), [Forceps](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
MATERIAL	Tungsten Carbide(TC)
AUTHORS	ADSON
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
MATERIAL	Tungsten Carbide(TC)
AUTHORS	ADSON
SPECIFICATIONS	serrated , straight