

TC SEMKEN TISSUE FORCEPS

straight; serrated; 1X2 teeth

18.0cm (7")

SKU: 021-240-18
Categories: [Open Surgery](#), [Forceps](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
MATERIAL	Tungsten Carbide(TC)
KLS Martin	12-240-18
AUTHORS	ADSON
SPECIFICATIONS	serrated , straight , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
MATERIAL	Tungsten Carbide(TC)
KLS Martin	12-240-18
AUTHORS	ADSON

SPECIFICATIONS

[serrated](#), [straight](#), [toothed](#)