

TC SEMKEN TISSUE FORCEPS

straight; serrated; 1X2 teeth

15.0cm (6")

SKU: 021-240-15
Categories: [Open Surgery](#), [Forceps](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
MATERIAL	Tungsten Carbide(TC)
KLS Martin	12-240-16
AUTHORS	SEMKEN
SPECIFICATIONS	serrated , straight , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
MATERIAL	Tungsten Carbide(TC)
KLS Martin	12-240-16
AUTHORS	SEMKEN

SPECIFICATIONS

[serrated](#), [straight](#), [toothed](#)