

TC SWEDISH MODELL NEEDLE HOLDER

straight; serrated

15.0cm (6")



SKU: 051-136-15
Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-652-15
AUTHORS	SCHWEDISCHES MODELL
Medicon	10.39.65
SPECIFICATIONS	serrated , straight

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-652-15

AUTHORS

[SCHWEDISCHES MODELL](#)

Medicon

10.39.65

SPECIFICATIONS

[serrated](#), [straight](#)