


TC VASCULAR NEEDLE HOLDER

straight; serrated

23.0cm (9")



SKU: 051-090-23/LH

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
AUTHORS	VASCULAR
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
AUTHORS	VASCULAR
SPECIFICATIONS	serrated , straight