

TC WANGENSTEEN DRESSING FORCEPS

straight; cross serrated

15.0cm (6")

SKU: 021-140-15
Categories: [Open Surgery](#), [Forceps](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD171
MATERIAL	Tungsten Carbide(TC)
KLS Martin	12-266-15
AUTHORS	WANGENSTEEN
Medicon	07.16.65
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD171

MATERIAL	Tungsten Carbide(TC)
KLS Martin	12-266-15
AUTHORS	WANGENSTEEN
Medicon	07.16.65
SPECIFICATIONS	serrated , straight