

## TC WANGENSTEEN NEEDLE HOLDER

straight; serrated

27.0cm (10.8")

**SKU:** 051-160-27

**Categories:** [Open Surgery](#), [Suture](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>Aesculap</b>	BM042
<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>KLS Martin</b>	20-656-27
<b>AUTHORS</b>	<a href="#">WANGENSTEEN</a>
<b>Medicon</b>	10.40.78
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>Aesculap</b>	BM042

<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>KLS Martin</b>	20-656-27
<b>AUTHORS</b>	<a href="#">WANGENSTEEN</a>
<b>Medicon</b>	10.40.78
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>