

TC WANGENSTEEN NEEDLE HOLDER

straight; serrated

27.0cm (10.8")

SKU: 051-160-27

Categories: [Open Surgery](#), [Suture](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM042
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-656-27
AUTHORS	WANGENSTEEN
Medicon	10.40.78
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM042

MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-656-27
AUTHORS	WANGENSTEEN
Medicon	10.40.78
SPECIFICATIONS	serrated , straight