

TC WEBSTER NADELHALTER

straight; smooth

13.0cm (5.2")




SKU: 051-029-13

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
AUTHORS	WEBSTER
SPECIFICATIONS	straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
AUTHORS	WEBSTER
SPECIFICATIONS	straight