

TC WEBSTER NEEDLE HOLDER

straight; serrated

13.0cm (5.2")



SKU: 051-026-13

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-604-13
AUTHORS	WEBSTER
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-604-13
AUTHORS	WEBSTER

SPECIFICATIONS

[serrated](#), [straight](#)