

TC WIRE CUTTING PLIER; ANGLED

soft wire up to Ø2.8mm; hard to Ø2.2mm

double convertible; 23.0cm (9.2")

SKU: 190-773-22

Categories: [Open Surgery](#), [Orthopaedics](#)



ZUSÄTZLICHE INFORMATIONEN


UNIT [pcs.](#)

TYPE [Wire Forceps](#)

Aesculap LX157

MATERIAL [Tungsten Carbide\(TC\)](#)

KLS Martin 22-523-22



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Wire Forceps](#)

Aesculap LX157

MATERIAL [Tungsten Carbide\(TC\)](#)

KLS Martin 22-523-22