

TC WIRE CUTTING PLIER; ANGLED

soft wire up to Ø2.8mm; hard to Ø2.2mm

double convertible; 23.0cm (9.2")


SKU: 190-773-22

Categories: [Open Surgery](#), [Orthopaedics](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Wire Forceps
Aesculap	LX157
MATERIAL	Tungsten Carbide(TC)
KLS Martin	22-523-22



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Wire Forceps
Aesculap	LX157
MATERIAL	Tungsten Carbide(TC)
KLS Martin	22-523-22